NEGATIVE PRESSURE WOUND THERAPY


Avance® Solo NPWT System

Manufacturer: Mölnlycke Health Care AB, Gamlestadsvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden

www.molnlycke.com/symbols

- STERILE EO**: The device is sterilized using ethylene oxide
- Do not reuse**
- Caution, see instructions for use**
- Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.**
- Follow instructions for use**
- Temperature limitation**
- Keep dry**
Keep away from rain
- Keep away from sunlight**
Keep away from heat
- REF**: Catalogue number
- Use by date / Expiry date**
- LOT**: Batch code
- SN**: Serial number
- MD**: Medical Device
- Do not use if package is damaged**
- Manufacturer**
- MR**: MRI unsafe
- Humidity limitation**
- Atmospheric pressure limitation**
- Applied part type BF**
- IP22**: Ingress Protection
- ETL Listed Marking**
- Separate collection of Waste Electrical and Electronic Equipment (WEEE)**
- System lasts up to 14 days**
- Leakage**
- Blockage**
- Low battery**

7.2. Avance Solo NPWT System – Normal use

Avance Solo Pump displays the following visual and audible signals to inform the user that the Avance Solo NPWT system is performing under normal operation.

AUDIBLE AND VISUAL DISPLAY	DESCRIPTION	OPERATION	COMMENT
	1. The pump activates for a short period. 2. All light indicators on the pump flash in turn. 3. The pump emits audible notifications: medium frequency followed by high frequency	Automatic self-check	Automatic self-check is performed when the batteries are correctly inserted in the pump and confirms that the pump is ready for use.
	1. At start-up of the pump, the green start button flashes once every second for 15 minutes. 2. During normal operation the green start button flashes two times every minute.	Therapy mode	1. To confirm that correct negative pressure is reached. 2. To confirm that the pump is operating correctly, and that therapy is maintained.
	1. The pump emits two short audible notifications 2. The audible notifications are repeated every 15 minutes as long as the therapy is paused.	Pause mode	To confirm that the pump and therapy have been paused. After 60 minutes the pump automatically restarts therapy.
	1. All light indicators flash with high intensity. 2. The pump emits three audible notifications: one high, one medium followed by one low frequency tone.	End of therapy	The therapy time of 14 days is complete.
	The pump emits short audible notification	Invalid button press	In case of an invalid button press

7.3. Avance Solo NPWT System – Alarms and trouble shooting

Avance Solo Pump displays the following visual signals and audible alarms to inform the user if risk for loss of therapy.

AUDIBLE AND VISUAL DISPLAY	DESCRIPTION	POSSIBLE CAUSE	COMMENTS TROUBLE SHOOTING
	1. The indication light for leakage flashes once per second. 2. If negative pressure is not maintained: • The light indicator for air leakage stays active, flashing once every second. • The pump repeatedly emits an audible alarm	Leakage Alarm Negative pressure is not being established due to an air leak in the system	To correct a leak, perform one or more of the following actions: Press around the dressing border to improve contact with the skin or if required add additional fixation strips around dressing borders. Ensure canister is securely attached to the pump. Ensure tubing is securely attached to the canister. Ensure that the dressing tubing is securely connected to the canister tubing.
	1. The indication light for blockage flashes once every second. 2. If the blockage condition persists • the light indicator for blockage stays active, flashing once every second	Blockage Alarm Negative pressure is not being established due to blockage	To correct a blockage, perform one or more of the following actions: Ensure that the tubing is not clamped. Ensure that tubing is not kinked. If the canister is full perform a canister change according to instructions for canister change.
	1. The indication light for battery flashes once every five (5) seconds when up to 24 hours of battery time remain. 2. When less than 4 hours of battery time remain • the light indicator for battery stays active, flashing once every second • the pump repeatedly emits an audible alarm	Low Battery Alarm	To correct a low battery alarm: Change the batteries according to instructions for battery change. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
	1. The indication lights for leakage, battery and blockage flash simultaneously once every second. 2. The pump repeatedly emits an audible alarm.	Internal Failure Alarm	The pump has an internal fault and cannot be started. Contact your healthcare professional or Mölnlycke Health Care.

8. Specifications Avance Solo Pump	
Nominal negative pressure	-125 mmHg
Maximum negative pressure	-150 mmHg
Mode of Operation	Continuous
Dimensions	Avance Solo Pump and Canister 50 ml 125x68x30 mm
Weight	Avance Solo Pump and Canister 50 ml < 130 g
Applied Part	Dressing, type BF
Battery	2x AA 1.5V Energizer L91
IP22	Ingress protection effective against fingers and similar objects. Protected against dripping water when tilted at 15°. Classification only valid when battery lid is closed.
Storage	Temperature 5°C/41 F to 25°C/77F, ambient humidity 10% to 75% non-condensing, ambient pressure 700 hPa to 1060 hPa
Transport	Temperature -35°C/-31 F to 63°C/145 F, ambient humidity 10% to 90% non-condensing, ambient pressure 700 hPa to 1060 hPa
Operation	Temperature 5°C/41 F to 40°C/104 F, ambient humidity 15% to 90%, non-condensing, ambient pressure 700 hPa to 1060 hPa
Low priority alarm signal, Alarm Volume 60 dBA	Leakage Alarm, Blockage Alarm, Low Battery Alarm, Internal Failure Alarm.
Information signals with lower priority than alarm signals	Pause mode, Therapy mode, Invalid button press, Pump Self-check, End of therapy, Leakage, Blockage, Low Battery.
Essential Performance	Activation of low priority alarms within two hours if degradation of Nominal negative pressure. Negative pressure not exceeding Maximum negative pressure longer than five minutes.

9. Safety
Avance Solo NPWT System complies with the General Requirements for Safety of Medical Electrical Equipment (IEC 60601-1). Avance Solo NPWT System is intended for home care use (IEC 60601-1-11).

10. Electromagnetic compatibility
Avance Solo Pump is tested in accordance with the requirements of IEC 60601-1-2. Exceeding test levels can cause failure to maintain Essential Performance.
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to Avance Solo Pump. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic Emissions		
Avance Solo Pump is intended for use in the electromagnetic environment specified below.		
EMISSION TESTS	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	Avance Solo Pump uses RF energy only for its internal function
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	Battery operated device
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Rx only Caution: Federal [US] Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Avance Solo NPWT System is intended to be used by healthcare professionals by directions of these instructions for use. Information for patient or lay person is provided in a separate patient user manual provided by Mölnlycke Health Care. The healthcare professional shall ensure that the patient user manual is handed to the patient or lay person as appropriate.

For prescription of therapy to home care, the prescribing healthcare professional should confirm that the patient or lay caregiver understands how the pump and canister works and how to operate on a daily basis. The prescribing healthcare professional should ensure that the patient or lay caregiver is able to perceive audible and visual notifications and alarms from the pump and to trouble shoot by guidance in Patient User Manual. The patient or lay person should be advised to contact prescribing healthcare professional if any concern on safe use of Avance Solo NPWT System.

- Product description**
Avance Solo Negative Pressure Wound Therapy (NPWT) System consists of Avance Solo Pump, Avance Solo Canister 50 ml, Avance Solo Border Dressing and Avance Solo Foam which together form a system for wound management via the application of negative pressure.
Avance Solo Pump, a battery powered single patient use pump with a 14 day lifespan, single button operated with visual and audible alarms and notifications
Avance Solo Canister 50 ml, a single use canister attached to the pump for collection of wound fluid and exudate
Avance Solo Border Dressing, a single use breathable soft silicone absorbent dressing with acrylic fixation strips
Avance Solo Foam, a single use polyurethane foam wound filler for cavity wounds
Avance Solo NPWT System maintains negative pressure nominally at -125 mmHg to the wound and enables exudate management by absorption and evaporation in Avance Solo Border Dressing. Excess exudate is collected in Avance Solo Canister 50 ml.
Avance Solo NPWT System is applicable on wound sizes (surface area x depth) up to 400 cm²/24 in² on low and moderately exuding wounds.
Avance Solo NPWT System is intended for use by healthcare professionals for therapy on patients in healthcare facilities and home care settings.
Avance Solo NPWT System is intended for adults.
- Indications for use**
Avance Solo NPWT System is indicated for patients who would benefit from wound management via the application of negative pressure, particularly as the device may promote wound healing through the removal of exudate, infectious material.
Avance Solo NPWT System is indicated for removal of low to moderate amounts of exudates from chronic, acute, traumatic, subacute and dehisced wounds, ulcers (such as diabetic, venous or pressure), surgically closed incisions, flaps, and grafts.
- Contraindications**
Avance Solo NPWT System is contraindicated for patients with the following conditions:
 - Malignancy in the wound or margins of the wound
 - Untreated and previously confirmed osteomyelitis
 - Non-enteric and unexplored fistulas
 - Necrotic tissue with eschar present
 - Exposed nerves, arteries, veins or organs
 - Exposed anastomotic site
- Warnings**
 - Application of NPWT may increase the risk of bleeding. If sudden or increased bleeding is observed, immediately take appropriate action to stop bleeding and seek emergency medical attention.
 - Patients at high risk of bleeding, such as patients receiving anticoagulant therapy or with altered haemostasis, must be monitored closely during therapy. Ensure to establish baseline before applying therapy.
 - Patients at risk of bleeding complications, due to e.g. a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels, should be monitored carefully during therapy.
 - Patients undergoing NPWT need frequent supervision. Routinely check that the negative pressure therapy is active, the dressing should be contracted and firm to the touch. If it is necessary to discontinue therapy, the time elapsed without negative pressure must comply with instructions from healthcare professional.
 - During therapy, make sure that the pump and tubings from the dressing and canister are positioned to eliminate the risk of
 - imprints
 - contamination
 - entrapment or strangulation
 - kinking or blocking of tubing
 - being exposed to sources of heat
 - Cover or remove sharp edges or bones in the wound prior to the application of the dressing due to the risk of puncturing of organs and blood vessels.
 - If defaturation is required, disconnect the pump and remove the dressing if dressing positioning interferes.
 - Avance Solo Pump is Magnetic Resonance (MR) unsafe, do not take into a Magnetic Resonance Imaging (MRI) environment. Avance Solo Border Dressings and Avance Solo Foam are MR safe. The impact of the dressing and foam on Magnetic Resonance Tomography (MRT)/Magnetic Resonance Imaging (MRI) imaging artifacts is unknown.
 - Consider use of non-adherent wound contact layer (Mepitel) to protect fragile tissue.
 - Do not place the pump with canister in water or other liquids. Disconnect the pump and canister if the ingress of water is observed.

- Avance Solo Pump shall not be used in oxygen rich environments where there is a danger of explosion e.g. a hyperbaric oxygen unit or in therapy/investigations involving microwaves. For the dressing, static electricity may be created when removing the release liner from the strips. Consider if the dressing needs to be removed for safety reasons.
- Avance Solo Pump is not suitable for use in the presence of flammable anesthetics.
- If the patient experiences sudden elevation in blood pressure or heart rate in response to stimulation of the sympathetic nervous system, immediately stop the therapy to help minimize sensory stimulation and seek emergency medical attention.
- To minimize the risk of bradycardia, do not place the dressing close to the vagus nerve.
- Do not use oxidizing agents such as hypochlorite solutions or hydrogen peroxide prior to use of the Avance Solo Foam.
- Products in Avance Solo NPWT System contain small parts which could be a potential choking hazard. Keep this device out of the reach of children.
- Keep this device out of the reach of pets.
- If canister or pump is broken, change the pump, replace with new product and restart therapy.
- 5. Precautions**
 - Before starting therapy, evaluate the patient's nutritional status and address severe malnutrition.
 - Signs of possible infection or complications must be addressed immediately. Monitor the device, wound, surrounding skin, and patient status accordingly to ensure effective and safe treatment and patient comfort.
 - For patients with ischemic condition or at application of a circumferential dressing, extra monitoring of wound status is required to avoid risk of compromised circulation.
 - Avance Solo Pump is equipped with visual and audible notifications and alarms. Ensure that carrying or placement of the pump allows the user to detect audible and visual notifications and alarms.
 - When Avance Solo Pump battery low alarm is triggered, replace the batteries in the pump. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
 - Fixation strips provided with the dressing shall only be applied on dressing borders. Do not apply fixation strips, or other occlusive dressings, across the dressing wound pad due to the risk of low breathability leading to maceration.
 - Application of certain skin protection products or using cleansing products prior to the application of the dressing, can affect the ability of the dressing and fixation strips to adhere securely and create sufficient sealing.
 - Do not use the products on patient and/or user with known hypersensitivity to the ingoing materials/components of the products.
 - Ingrowth of tissue may occur if the dressing or wound filler is not changed in accordance with recommendations or as appropriate for the wound condition of the individual patient (See Section 6.6).
 - Wounds deeper than 0.5 cm (0.2 inches) are likely to require Avance Solo Foam as wound filler, to ensure dressing to establish intimate contact with all areas of the wound bed.
 - For dressing applications with the larger sizes of dressing, the wound should not generally be deeper than 2 cm (0.8 inches). If excessive exudate is present, consider the use of traditional NPWT until the level of exudate is low to moderate and then transition the patient to Avance Solo NPWT System.
 - Avance Solo NPWT System may be applied in conjunction with compression therapy. Make sure not to cover the transfer port on the dressing to reduce to risk of imprints.

- 6.2. Pump Set up**
 1. Attach the canister to the pump by pushing the canister until it clicks on both sides.
 2. Fit the canister tubing into the holder on the back of the pump.
 3. Insert the batteries into the pump battery compartment. Ensure that the positive terminal (marked +) and negative terminal (marked -) of each battery matches the +/- label in the battery compartment. Close the battery compartment by sliding the lid back in place.
 4. When the batteries are correctly inserted, the pump performs automatic self-check.
 5. After completed self-check, the pump will remain paused until started. Audible notification is repeated every 15 minutes as long as the therapy is paused. To start therapy press and hold the green pump start button, release after two (2) seconds.

If not started directly, the pump will start automatically after 60 min.
When batteries are inserted into the pump for the first time, an internal timer starts. Document date and time for start of therapy in the patient's notes. Avance Solo Pump has a lifetime of 14 days from the first insertion of the batteries. It is not possible to restart the pump after therapy time has ended.
- 6.3. Dressing Application Procedure**
To apply, use clean/aseptic or sterile techniques in accordance with local protocol.
1. Cleanse and debride the wound bed as instructed by healthcare professional.
2. Cleanse the peri wound skin and pat dry.
3. Eliminate or cover sharp edges and bone fragments with non-adherent wound contact layer due to risk of puncturing organs or blood vessels while under negative pressure. Document the use of wound contact layer in the patient's notes.
4. If a wound filler is used, see Section 6.5 Application with wound filler.
5. Decide in which direction the dressing should be placed to avoid bends and kinks of the tubing and to avoid discomfort to the patient. Depending on the patient's primary position the transfer port on the dressing shall be placed upmost from the wound.
6. Grip the center part of the release film of the dressing and pull to expose the adhesive surface.
7. Without stretching, position the dressing centrally over the wound and ensure that dressing borders are placed on intact skin.
8. Gently remove the remaining release films of the dressing. Start with the longer portion of the release film, away from the tubing. Do not stretch the product during application. Smooth down the dressing to remove any creases and press gently to ensure contact between the dressing and wound bed.
9. Connect the canister tubing to the dressing tubing using the connectors on the end of each tubing. Ensure that the tubings are not clamped.
10. To start therapy, press and hold the green pump start button, release after two (2) seconds so that the pump activates and the green start button flashes. Negative pressure shall be achieved within two (2) minutes after pump start.

- 6.7. Dressing removal procedure**
To remove the dressing, perform the following steps:
27. If the pump is active, pause therapy by pressing down the green start button, release after two (2) seconds.
28. Block both the canister tubing and dressing tubing by positioning the slide clamps next to the connector and slide across the tubings until secured.
29. Disconnect the canister tubing from the dressing tubing by squeezing the connector from both sides and pull apart.
30. Gently peel back one corner of the film of fixation strips and dressing and stretch to facilitate breakage of the seal. Proceed with this technique (in the direction of hair growth) until the film is completely removed.
31. Remove the dressing by pulling it in the direction of the wound, not across the wound. Pull off the dressing following the surface all the way, without lifting the dressing perpendicular to the wound.
32. If wound filler is used, remove the wound filler gently. If adherence of wound filler to wound bed is observed, consider moistening the filler material. Make sure no fragments are left in the wound when the dressing is changed.
If the patient experiences pain during dressing removal, consider the use of pain relief medication, as instructed by healthcare professional, when changing dressings. Consult patient's notes to ensure that all materials used have been removed.
- 6.8. Changing the Canister**
The need to change canister is detected either by visual inspection of canister fill level through the transparent window on the back of the canister or by pump blockage alarm.
To replace the canister, perform the following steps:
33. If the pump is active, pause therapy by pressing down the green start button, release after (2) seconds.
34. Clamp the canister tubing and dressing tubing by positioning the slide clamps next to the connector and slide across the tubings until secured. Blocking the tubings minimizes fluid leakage when disconnecting the dressing from the canister.
35. Disconnect the canister tubing from the dressing tubing by squeezing the connector from both sides and pull apart.
36. Remove the canister by pushing the spring buttons on both sides and pull.
37. Attach a new canister to the pump by pushing the canister until it clicks on both sides. To continue therapy, connect the canister tubing to the dressing tubing, release the clamp on the dressing tubing and restart the pump by pressing down the green start button, release after two (2) seconds.
- 6.9. Changing the Batteries**
The pump indicates when the batteries are low, as described in Section 7. Avance Solo NPWT System Indicators, Alarms and Trouble shooting. Batteries shall be replaced once pump alarms for Low battery or after 7 days. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
To replace the batteries, perform the following steps:
38. If pump is still active, pause the pump by pressing down the green start button, release after two (2) seconds.
39. Open the battery compartment by sliding the lid. Remove the batteries. Insert new batteries, ensuring that the positive terminal (marked +) and negative terminal (marked -) of each battery matches the +/- label in the battery compartment. Close the battery compartment lid.

Electromagnetic Immunity			
Avance Solo Pump is intended for use in the electromagnetic environment specified below.			
IMMUNITY TESTS	BASIC EMC STANDARD OR TEST METHOD	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	Home healthcare environment
Electrostatic Discharge	61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Radiated RF EM fields	61000-4-3		10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment	61000-4-3	30 cm minimum separation distance from radio transmitter	
Rated power frequency magnetic fields	61000-4-8	30 A/m 50 Hz or 60 Hz	


- 11. Caution**
Avance Solo NPWT System must be used in accordance with these Instructions for Use. Read these instructions before using the system and have them available during use. Failure to read and understand these instructions may lead to misuse of the system and improper performance. These instructions are a general guide for the use of the product. Specific medical situations must be addressed by a clinician.
- 12. Other information**
If any serious incident has occurred in relation to Avance Solo NPWT System, it should be reported to Mölnlycke Health Care.

7.1. Avance Solo NPWT System Indicators	
DISPLAY	LIGHT INDICATOR
	Leakage
	Blockage
	Low Battery

Australian sponsor address : Mölnlycke Health Care Pty Ltd, 12 Narabang Way, Belrose, NSW 2085, AUSTRALIA

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