EN **Instructions for Use PPE Products for Biogel Surgeons, Biogel Super-Sensitive, Biogel Indicator System, Biogel Eclipse, Biogel Eclipse Indicator System, Biogel Eclipse Indicator Underglove, Biogel Skinsense, Biogel Skinsense Indicator System, Biogel Skinsense Indicator Underglove, Biogel PI UltraTouch, Biogel PI Indicator System, Biogel PI Indicator Underglove, Biogel PI Micro, Biogel PI Micro Indicator System, Biogel PI Micro Indicator Underglove, Biogel PI UltraTouch S, Biogel PI UltraTouch S Indicator System, Biogel PI UltraTouch S Indicator Underglove.**

These gloves are in conformity with Regulation (EU) 2016/425 (Personal Protective Equipment). The gloves comply with EN420:2003+A1:2009.

USE:
For Single Use Only
Regulation (EU) 2016/425 category III product of complex design. The gloves are designed to provide a barrier to microorganisms. The gloves are liquid resistant and can be used for splash protection against certain chemicals. Please consult Mölnlycke Health Care if more information is required on the level of performance of the gloves with specific chemicals. Chemical data is available from Mölnlycke Health Care.

The actual duration of protection provided in the workplace may vary considerably from these performance levels due to other factors influencing the performance, such as temperature, abrasion, degradation.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defect or imperfections.

The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PRECAUTIONS:
Gloves damaged in use should be removed and replaced as soon as damage is detected.

DISPOSAL:
Used gloves/sleeves may be contaminated with infectious or other hazardous materials. Dispose of according to Local Authority Regulations.

EU type examination certificate from BSI Group (No 2797) and is the notified body for Regulation (EU) 2016/425: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

For Declaration of Conformity, please refer to www.molnlycke.com/Declaration-of-Conformities

EN ISO 374-5:2016	EN ISO 374-1:2016 Type C
	
VIRUS AQL 0.65	K

The gloves are available in size 5.5 to 9 (Biogel Surgeons, Biogel Super-Sensitive, Biogel Eclipse, Biogel Eclipse Indicator Underglove, Biogel Skinsense, Biogel Skinsense Indicator Underglove, Biogel PI UltraTouch, Biogel PI Indicator Underglove, Biogel PI Micro, Biogel PI Micro Indicator Underglove, Biogel PI UltraTouch S, Biogel PI UltraTouch S Indicator Underglove)

The gloves are available in size 5.5 to 8.5 (Biogel Eclipse Indicator System, Biogel PI Micro Indicator System, Biogel Indicator System, Biogel PI Indicator System, Biogel Skinsense Indicator System, Biogel PI UltraTouch S Indicator System)

EN **Gebrauchsanweisung für PSA-Produkte für Biogel Surgeons, Biogel Super-Sensitive, Biogel Indicator System, Biogel Eclipse, Biogel Eclipse Indicator System, Biogel Eclipse Indicator Underglove, Biogel Skinsense, Biogel Skinsense Indicator System, Biogel Skinsense Indicator Underglove, Biogel PI UltraTouch, Biogel PI Indicator System, Biogel PI Indicator Underglove, Biogel PI Micro, Biogel PI Micro Indicator System, Biogel PI Micro Indicator Underglove, Biogel PI UltraTouch S, Biogel PI UltraTouch S Indicator System, Biogel PI UltraTouch S Indicator Underglove.**

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