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
Referanser:
 1. Upton D, Solowiej K. Pain and stress as contributors to delayed wound healing. *Wound Practice and Research* 2010, Vol. 18(3).
 2. Gee Kee EL, Kimble RM, Cuttle L, Khan A, Stockton KA. Randomized controlled trial of three burns dressings for partial thickness burns in children. *Burns* 2015, 41(5):946-955.
 3. Silverstein P, Heimbach D, Meites H et al. An open, parallel, randomized, comparative, multicenter study to evaluate the cost-effectiveness, performance, tolerance, and safety of a silver-containing soft silicone foam dressing (intervention) vs silver sulfadiazine cream. *J Burn Care Res* 2011, 32(6): 617-626.
 4. Gee Kee EL, Stockton K, Kimble RM et al. Cost-effectiveness of silver dressings for paediatric partial thickness burns: An economic evaluation from a randomized controlled trial. *Burns* 2017, 43(4): 724-732.
 5. Agarwal S, Harish V, Roberts S et al. Treatment of partial thickness burns: a prospective, randomised controlled trial comparing Biobrane, Acticoat, Mepilex Ag and Aquacel Ag. *J Burn Care Res* 2020, 42(5): 934-43.
 6. Tang H, Lv G, Fu J et al. An open, parallel, randomized, comparative, multicenter investigation evaluating the efficacy and tolerability of Mepilex Ag versus silver sulfadiazine in the treatment of deep partial-thickness burn injuries. *J Trauma Acute Care Surg* 2015, 78(5): 1000-1007.
 7. Gleeson AP, Stuart M, Wilson B, Phillips B. Ultrasound assessment and conservative management of inversion injuries of the ankle in children: Plaster of Paris versus tubigrip. *Journal of Bone and Joint Surgery - Series B* 1996, 78(3):484-7.
 8. Naem M, Rahimniaj M, Rahimniaj NA, Idrees Z, Shah GA, Abbas G. Assessment of functional treatment versus plaster of Paris in the treatment of grade 1 and 2 lateral ankle sprains. *Journal of Orthopaedics and Traumatology* 2014, 16(1):41-6.
 9. Mölnlycke Health Care. Data on file. Tubigrip. 10. Eytier C, Gazeau E, Beneteau G, Verfaillie G. Convenience and tolerance of the combination of a soft silicone foam dressing and a two-way stretch tubular bandage in the management of local wounds. *Journal des plaies et cicatrisations* 2013, 18(8):38-44.
 11. Mölnlycke Health Care. Data on file. Tubifast. 12. Mölnlycke Health Care. Data on file. Mepiform. 13. Wigger-Albert W, Kuhlmann M, Wilhelm D, Mrowietz U, Eichhorn K, Ortega J, et al. Efficacy of a polyurethane dressing versus a soft silicone sheet on hypertrophic scars. *Journal of wound care* 2009, 18(5):208, 10-4.
 14. Armstrong DG, Bohn G, Glaz P, Kavros SJ, Kriner R, Snyder R et al. Expert recommendations for the use of hypochlorous solution: science and clinical application. *Dermatol Wound Manage* 2015, 61(5 suppl): 45-185.
 15. Consensus on Wound Antisepsis. Update 2018. *Skin Pharmacol Physiol* 2018;31:28-58.
 16. Chadwick P, Taherinejad F, Hamburg K, Waring M. Clinical and scientific data on a silver-containing soft-silicone foam dressing: an overview. *Journal of Wound Care* 2009, 18(11):483-491.
 17. Mölnlycke Health Care. Data on file. Mepilex Transfer Ag. 18. Schweiger H, Smith D, Cruse CW et al. An open, non-controlled, single-centre, clinical investigation to evaluate efficacy when using a soft silicone wound contact layer containing silver. Poster presentation at the 15th European Burns Association Congress, Vienna, Austria, 2013.
 19. Mölnlycke Health Care. Data on file. Mepitel Ag. 20. Mölnlycke Health Care. Data on file. Melgisorb Plus. 21. Chadwick P, McCauley J. Open, non-comparative, multicenter post-market clinical study of the performance and safety of a gel-like fibre wound dressing on diabetic foot ulcers. *Journal of Wound Care* 2016, 25(4): 290-300.
 22. Smet S, Beebe H, Saine L, Suis E, Henricks B. Open, non-comparative, multi-centre post-market clinical follow-up investigation to evaluate performance and safety on pressure ulcers when using a gel-like fibre dressing as intended. Poster Presentation at European Pressure Ulcer Advisory Panel Conference 2015, Ghent, Belgium.
 23. Mölnlycke Health Care. Data on file. Exufiber. 24. David F et al. A randomised, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidcolloid wound contact layer (UrgoTul) in the treatment of acute wounds. *International Wound Journal* 2018, 25. Mölnlycke Health Care. Data on file. Avance Solo. 26. Nischwitz SP, Luze H, Popp D, Winter R, Draschl A, Schellnegger M, Kargl L, Rappi T, Giretzlehner M, Kamolz LP. Global burn care and the ideal burn dressing reloaded—A survey of global experts. *Burns* 2021, 47: 1665-1674.

Les mer på www.molnlycke.no
 Mölnlycke Health Care AS, postboks 6229, Etterstad, 0603 Oslo. Tlf: 22 70 63 70.
 The Mölnlycke, Mepilex, Mextra, Tubifast, Exufiber, Mepitel, Mepiform, Mesorb, Granudacyn, Tubigrip, Melgisorb, Avance and Safetac trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2023 Mölnlycke Health Care AB. All rights reserved. 2023NOWC3742305

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