Venøse leggsår

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Referanser: 1. Strohal R., Gerber V., Kröger K., et al. Expert consensus on practical aspects of wound therapy with hemoglobin spray. *Wound management* 2016; 5: 276-284. 2. Bjørk, R., Ehmann S., S.T.R.I.D.E. Professional guide to compression garment selection for the lower extremity. *Journal of Wound Care* 2019; 28(16 suppl.1):1-44. Consensus recommendations. *Wounds International* 2015; 3. Franks, P., Barker, J., Collier, M., et al. Management of patients with venous leg ulcers: challenges and current best practice. *J Wound Care*, 25, 6, Suppl, 1-67, 2016. 4. Harding, K., Dowsett, C., Fias, L., et al. Simplifying venous leg ulcer management. 5. Teagle, A. and Hargest, R. Management of pyoderma gangrenosum. *Journal of the Royal Society of Medicine* 2014; 107(6):228-236. 6. White R. Evidence for atraumatic soft silicone wound dressing use. *Wounds UK* 2005;1(3):104-109. 7. Waring, M., Bielefeld, S., Matzold, K.P., Butcher, M. An evaluation of the skin stripping of wound dressing adhesives. *J Wound Care* 2011;20:412-22. 8. White, R. A multinational survey of the assessment of pain when removing dressings. *Wounds UK* 2008;4:14-22. 9. David, F., Wutze, J.-L., Breton, N., et al. A randomized, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidocolloid wound contact layer (UrgoTul) in the treatment of acute wounds. *Int Wound J* 2017 doi:10.1111/iwj.12853. 10. Gotschall, C.S., Morrison, M.I., Eichelberger, M.R. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. *J Burn Care Rehabil* 1998;19:279-83. 11. Silverstein, P., Heimbach, D., Meites, H., et al. An open, parallel, randomized, comparative, multicenter study to evaluate the cost effectiveness, performance, tolerance, and safety of a silver-containing soft silicone foam dressing (Intervention) vs silver sulfadiazine cream. *J Burn Care Res* 2011;32:617-26. 12. Gee Kee, E.L., Kimble, R.M., Cuttle, L., Khan, A., Stockton, K.A. Randomized controlled trial of three burn dressings for partial thickness burns in children. *Burns* 2015;41:946-55. 13. Meaume, S., Van De Loooverbosch, D., Heyman, H., Romanelli, M., Gangherotti, A., Charpin, S. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. *Ostomy Wound Manage* 2003; 49(9): 44-51. 14. Wiberg, A.-B., Feili, F., Daun, E.-K. Preventing maceration with a soft silicone dressing: *in vitro* evaluation. Poster presentation at the 3rd Congress of the World Union of Wound Healing Societies, Toronto, Canada, 2008. 15. Bredow J, Hoffmann K, Helmrich M, Eysel P, Zaroghooni K. Randomized clinical trial to evaluate performance of flexible self-adherent absorbent dressing coated with silicone layer after hip, knee or spinal surgery in comparison to standard wound dressing. Poster presentation at the 5th Congress of the World Union of Wound Healing Societies, Florence, Italy, 2016. 16. Mölnlycke Health Care data on file, report 20160105-002. 17. Lantini, A., Diegel, C., Scheske, J., Schmitt, C., Brönner, A., Jodl, H. Mepilex XT in practice: results of a study in German specialist wound care centres. *Wounds International* 2015; 6(4):18-22. 18. ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Knees. Data on file. 19. ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Elbows. Data on file. 20. ALTEN Finite Element Modelling simulation. Laboratory report no. PD-530246. 21. Mölnlycke Health Care. Mepilex® Border Flex-Bacteria encapsulation. Report no. PD-537072. Data on file. 22. Chadwick, P., Tahernejad, F., Hamburg, K., Waring, M. Clinical and scientific data on a silver-containing soft-silicone foam dressing: an overview. *J Wound Care* 2009;18:483-91. 23. External lab report: NAMS A 09C 29253 01/09C 29253 02. 24. Kies C.L., Murrah, C.P., Smith, K., et al. Achieving and sustaining zero. Preventing surgical site infections after isolated coronary artery bypass with saphenous vein harvest through implementation of a staff-driven quality improvement process. *Dimensions Crit Care Nurs* 2015;34:265-72. 25. Smet, S., Beele, H., Saine, L., Suijs, E., Hennickx, B., Open, noncomparative, multi-centre post market clinical follow-up investigation to evaluate performance and safety on pressure ulcers when using a gelling fibre dressing as intend-ed. Poster Presentation at European Pressure Ulcer Advisory Panel Conference, 2015, Ghent, Belgium. 26. Chadwick P, McCordle J, Open, non-comparative, multicentre post clinical study of the performance and safety of a gelling fibre wound dressing on diabetic foot ulcers. *J Wound Care* 2016;25:290-300. 27. Davies P, McCarty S. An in-use product evaluation of a gelling fibre dressing in wound management. E-poster presentation at Wounds UK Conference, 2017, Harrogate, United Kingdom. 28. Lev-Tov et al. An interim analysis of clinical investigation to evaluate exudate management and comfort of use of an antimicrobial gelling fibre dressing in medium to highly exudative wounds. Poster presented at the Symposium of Advanced Wound Care, Fall meeting 2018, Las Vegas, NV, USA. 29. Mölnlycke Health Care Laboratory Report PD-520425 (unpublished). Retention 30. Mölnlycke Health Care. Exufiber. Gesellschaft für Versorgungskonzepte in der Wundbehandlung (GWW) mbH, Stuttgart, Germany. Data on file (unpublished report), 2017. 31. Mölnlycke Health Care Laboratory Report PD-521248 (unpublished). Tensile strength. 32. Mölnlycke Health Care Laboratory Report PD-521232 (unpublished). Absorption under compression. 33. Mölnlycke Health Care Laboratory Report PD-522900 (unpublished). Absorption of thick exudate 34. Mölnlycke Health Care Laboratory Report PD-521245 (unpublished). Absorption of blood. 35. Hamburg K et al. Antimicrobial effect of a new silver-containing gelling fibre dressing against common wound pathogens. Poster presented at the Symposium on Advanced Wound Care Spring meeting/ Wound Healing Society (WHS) Annual Meeting 2017, Apr 05-09, 2017, San Diego, CA, USA. 36. Gil J et al. 2017. Evaluation of a Gelling Fiber Dressing with Silver to Eliminate MRSA Biofilm Infections and Enhance the Healing. Poster presented at the Symposium on Advanced Wound Care Spring meeting/ Wound Healing Society (WHS) Annual Meeting 2017, Apr 05-09, 2017, San Diego, CA, USA. 37. Valdes et al. 2017. Evaluation of a Gelling Fiber Dressing with Silver to Eliminate P. a. Biofilm Infections and Enhance the Healing. Poster presented at the Symposium on Advanced Wound Care Spring meeting/ Wound Healing Society (WHS) Annual Meeting 2017, Apr 05-09, 2017, San Diego, CA, USA. 38. External lab report: NAMS A 11C, 51788, 01. 39. Grocott Patricia. Clinical investigation of Mepilex® Transfer. Clinical investigation of a silicone dressing in product development phase in the palliative management of patients with pressure sores and malignant wounds. study id M1N101 UK, London UK, 2000. 40. Tickle, J., Fletcher, J. Mextra Superabsorbent made easy. *Wounds UK* 2012;8: 1-4. 41. Eming, S., et al. (2008) The inhibition of matrix metalloproteinase activity in chronic wounds by polyacrylate superabsorbent Biomaterials; 29(19):2932-40. 42. Wiegand, C., et al. Polyacrylate superabsorbents bind inflammatory protease *in vitro*, Poster at Wounds UK, 2008. 43. Wiegand C., et al. (2013) A superabsorbent polymer-containing wound dressing efficiently sequesters MMPs and inhibits collagenase activity *in vitro*. *Journal of Materials Science Materials in Medicine*. 2013;24(10):2473-2478. 44. Hunt, S., Elg, F. Clinical effectiveness of hemoglobin spray (Granulox®) as adjunctive therapy in the treatment of chronic diabetic foot ulcers. *Diabetic Foot & Ankle* 2016; 7: 33101. 45. Elg, F., Hunt, S. Hemoglobin spray as adjunctive therapy in complex wounds: meta-analysis versus standard care alone in pooled data by wound type across three retrospective cohort controlled evaluations. *SAGE Open Med*. 2018 Jun; 6:2050312118784313. 46. Hunt, S., Elg, F. The clinical effectiveness of hemoglobin spray as adjunctive therapy in the treatment of chronic wounds. *J Wound Care*. 2017 Sep; 26(9):558-568. 47. *In vitro* suspension test [EN13727, EN 13624, EN 13704, EN 14476 - phase 2] with Granudacyn® wound irrigation solution. 48. Consensus on Wound Antiseptics: Update 2018, *Skin Pharmacol Physiol* 2018;31:28-58. DOI: 10.1159/000481545.

Dette er ment som veiledning og kan ikke erstatte klinisk vurdering. Enhver kliniker er ansvarlig for omfattende vurdering og individuell behandlingsplan tilpasset pasientens behov.