


Mölnlycke® Pulsed Lavage

Truly customised orthopaedic solutions



Battery version



AC version

Mölnlycke Pulsed Lavage is DEHP free

Our single-use high pressure pulsed lavage system effectively cleans bone and removes debris during total joint procedures. This increases cement penetration into the cancellous bone for potentially stronger implant fixation and reduced revision rates. It is also an effective mechanical wound debridement system for cleaning wounds, removing dirt, necrotic tissue and debris.

Available in two versions

- **AC adapter:** reduced environmental impact thanks to lack of alkaline batteries. Single-use handpiece and re-usable AC adapter
- **Built-in battery:** a flexible and mobile version with battery pack integrated into the handpiece, for quick and easy battery removal

All included means no hidden costs

- Ergonomic and comfortable handpiece
- Short fan spray tip with double splash shield
- Long femoral canal tip
- Integrated suction and irrigation
- Cleaning rod for solving potential occlusion in the femoral tip

Single packed or included in Mölnlycke® Procedure trays

- Include Pulsed Lavage in your customised procedure-specific tray for efficiency gains
- For hip revision surgery we have a special tip to handle small bone and cement parts. Available as sterile single packed
- Also available as a sterile single packed

Filtered revision Tip

- For hip revision surgery we have a special tip to handle small bone and cement parts. Available as sterile single packed

* Pulsed Lavage's legal manufacturer is Kaiser Technology co., Ltd. Mölnlycke Health Care is importer and distributor.



Mölnlycke®

Mölnlycke® Pulsed Lavage – AC version

Sterile single packed

Ref. No.	Description
9750000	Pulsed Lavage System, AC Power
9752000	AC Power Supply (AC100~240V), Europe
9752010	AC Power Supply (AC100~240V), DK
9752020	AC Power Supply (AC100~240V), UK
9752030	AC Power Supply (AC100~240V), AUS
9752040	AC Power Supply (AC100~240V), CH
9751000	Coaxial Canal Brush Tip
9751001	Cone Tip
9751002	Coaxial Canal Tip
9751003	Hip Revision Tip

Available in trays

Ref. No.	Description
2322130-00	Pulsed Lavage System, AC Power
2320481-00	Coaxial Canal Brush Tip



AC Version




AC Adapter




Coaxial Canal Brush Tip



Hip Revision Tip



Cone Tip



Coaxial Canal Tip


Mölnlycke® Pulsed Lavage – Build-in Battery

Sterile single packed


Ref. No.	Description
9750003	Pulsed Lavage System, Built-in Battery
9751004	Coaxial Canal Brush Tip
9751005	Cone tip
9751006	Coaxial Canal Tip
9751007	Hip Revision Tip

Available in trays


Ref. No.	Description
2324050-00	Pulsed Lavage System, Built-in Battery
2325431-00	Coaxial Canal Brush Tip




Built-in Battery




Coaxial Canal Brush Tip



Cone Tip



Hip Revision Tip



Coaxial Canal Tip

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00.
The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2021 Mölnlycke Health Care AB. All rights reserved. HQIM002440

