

## PRODUCT DATA SHEET

### Tubifast with 2 way stretch technology 3 metre

Product name	Tubifast with 2 way stretch technology				
Product description	Elasticated viscose tubular bandage				
Intended use	For dressing retention, the treatment of atopic eczema as part of the wrapping technique and as an undercast stockinette				
Product Codes	Art.No	Product Name	Size	Pcs/Shelf cont	Pcs/TRP
	992485-02	Tubifast Green	3m	6	168
	992486-02	Tubifast Blue	3m	6	96
	992487-02	Tubifast Yellow	3m	6	96
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden				
Country of origin	UK				
Product classification according to MDD	Class I				
CE-marked	Yes, Notified Body BSI 0086				
Products material	Viscose and elastane covered with polyamide				
Durability	3 years				
Packaging	The purpose of the shelf container and transport box is to contain and protect material from dust, moisture and damage during transport, handling and storage.				
EAN Codes	Art No	TRP	RET	CON	
	992485-02	5055158009066	5055158009059	5055158009042	
	992486-02	5055158009103	5055158009097	5055158009080	
	992487-02	5055158009141	5055158009134	5055158009127	