

Tubigrip®

Multi-purpose tubular support bandage

Provides firm support for sprains, strains and swelling



What size to choose?

Product code	Tubigrip size	Minimum limb circumference (cm)	Minimum limb circumference (inch)	Maximum limb circumference (cm)	Maximum limb circumference (inch)	Sub-bandage pressure (mmHg)
1435	A	10	3.9	14.5	5.7	20 - 30
1436	B	14.5	5.7	29	11.4	10 - 20
1443	C	22	8.7	35	13.8	
1437	D	26	10.2	40	15.7	
1434	E	29	11.4	46	18.1	
1438	F	35	13.8	53	20.9	
1439	G	40	15.7	62	24.4	5 - 15
1440	J	54	21.3	82	32.3	

Tubigrip is available under a full range of sizes to suit most applications. The sub-bandage pressures are obtained in double layer.

Reference: Mölnlycke Health Care. Data on file.

Find out more at www.molnlycke.ca

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00. The Mölnlycke, Tubigrip trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2020 Mölnlycke Health Care AB. All rights reserved. HQIM002264.

