

MÖLNLYCKE - WOUND CARE PRODUCTS



PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

BORDERED FOAM DRESSINGS

	<p>Mepilex® Border Flex Absorbent self-adherent foam dressing with 360° stretch Flex Technology</p>  	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds • Diabetic foot ulcers • Leg and foot ulcers • Traumatic wounds, such as skin tears • Pressure injury prevention* and treatment <p>* When a repositioning/offloading protocol is used</p>	<p>7.5cm x 7.5cm 10cm x 10cm 12.5cm x 12.5cm 15cm x 15cm 15cm x 20cm</p>	<p>595211 595311 595011 595411 595611</p>
	<p>Mepilex® Border Flex Oval Absorbent self-adherent foam dressing with Flex Technology</p>  	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds for challenging body contours • Diabetic foot ulcers • Leg and foot ulcers • Traumatic wounds, such as skin tears • Pressure injury prevention* and treatment <p>* When a repositioning/offloading protocol is used</p>	<p>7.8cm x 10cm 13cm x 16cm 15cm x 19cm</p>	<p>583500 583300 583400</p>
	<p>Mepilex® Border Heel Absorbent self-adherent foam dressing for the heel area</p>   <p>Added benefits:</p> <ul style="list-style-type: none"> • Thicker borders to minimise rolling and to allow efficient handling • Optimised seal for enhanced protection of high-risk heel area. • Enhanced coverage of high risk areas • Additional handling tabs 	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds for the heel area • Pressure injury prevention and treatment • Diabetic foot ulcers • Heel ulcers • Traumatic wounds 	<p>22cm x 23cm</p>	<p>282750</p>
	<p>Mepilex® Border Sacrum Absorbent self-adherent foam dressing for the sacral area</p>   <p>Added benefits:</p> <ul style="list-style-type: none"> • Thicker borders to minimise rolling and to allow efficient handling • Optimised gluteal seal for enhanced protection of high-risk sacral area. • Enhanced coverage of high risk areas • Additional handling tabs 	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds for the sacral area • Sacral pressure injury prevention and treatment (all stages) • Protection from moisture, shear, friction 	<p>16cm x 20cm 22cm x 25cm</p>	<p>282050 282450</p>
	<p>Mepilex® Border Lite Thin absorbent self-adherent foam dressing</p>  	<ul style="list-style-type: none"> • Non to low exuding wounds • Foot ulcers and diabetic foot ulcers • Leg ulcers • Traumatic wounds • Skin tears • Wounds with compromised and / or fragile skin 	<p>4cm x 5cm 5cm x 12.5cm 7.5cm x 7.5cm 10cm x 10cm 15cm x 15cm</p>	<p>281000 281100 281200 281300 281500</p>

MÖLNLYCKE - WOUND CARE PRODUCTS












PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

NON-BORDERED FOAM DRESSINGS

	<p>Mepilex® XT Highly conformable absorbent foam dressing, with specially designed exudate channels to assist in the absorption of thick exudate</p>  	<ul style="list-style-type: none"> • Low to moderately exuding wounds in all healing stages • Chronic wounds with thick exudate • Leg ulcers and diabetic foot ulcers • Pressure injuries • Traumatic wounds • Can be used under compression 	<p>5cm x 5cm 10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>211015 211100 211200 211300 211400 211500</p>
	<p>Mepilex® Highly conformable and absorbent foam dressing</p>  	<ul style="list-style-type: none"> • Low to moderately exuding wounds • Leg ulcers and diabetic foot ulcers • Pressure injuries • Traumatic wounds • Partial thickness burns • Can be used under compression 	<p>5cm x 5cm 10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>294015 294100 294200 294300 294400 294500</p>
	<p>Mepilex® Heel Absorbent foam dressing for the heel area</p>  	<ul style="list-style-type: none"> • A wide range of exuding wounds on the heel or malleolus • Foot ulcers / diabetic foot ulcers 	<p>13cm x 20cm</p>	<p>288100</p>
	<p>Mepilex® Lite Thin absorbent foam dressing</p>  	<ul style="list-style-type: none"> • Non to low exuding wounds • Diabetic foot ulcers • Leg ulcers • Radiation skin reactions • Areas requiring protection of compromised and / or fragile skin 	<p>6cm x 8.5cm 10cm x 10cm 15cm x 15cm 20cm x 50cm</p>	<p>284000 284100 284300 284500</p>
	<p>Lyof foam® Max Polyurethane foam dressing</p> 	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds • Leg ulcers • Pressure injuries 	<p>7.5cm x 8.5cm 10cm x 10cm 15cm x 15cm 20cm x 20cm</p>	<p>603200 603201 603204 603206</p>

WOUND CONTACT LAYERS

	<p>Mepitel® One Wound contact layer dressing with SafetaC® on one side</p>  	<ul style="list-style-type: none"> • Skin tears and skin abrasions • Lacerations • Surgical incisions • Fixation of partial and full thickness grafts • Partial thickness burns • Diabetic foot ulcers • Venous and arterial leg ulcers <p>* Sizes with pre-printed symbols</p>	<p>5cm x 7.5cm 6cm x 7cm 9cm x 10cm 7.5cm x 10cm 13cm x 15cm 10cm x 18cm 17cm x 25cm 24cm x 27.5cm 27.5cm x 50cm 9.5cm x 1.5cm</p>	<p>289100 289170* 289270* 289300 289470* 289500 289700 289670 289750 289800</p>
	<p>Mepitel® Wound contact layer dressing with SafetaC® on both sides</p>  	<ul style="list-style-type: none"> • Skin tears and skin abrasions, lacerations • Surgical incisions • Partial thickness burns • Fixation of partial and full thickness grafts • Diabetic foot ulcers, venous and arterial leg ulcers 	<p>5cm x 7.5cm 7.5cm x 10cm 10cm x 18cm 20cm x 30cm</p>	<p>290510 290710 291010 292005</p>
	<p>Mepilex® Transfer Exudate transfer foam dressing</p>  	<ul style="list-style-type: none"> • Highly exuding wounds, venous/arterial leg ulcers • Difficult-to-dress areas • Fungating tumours • Donor sites 	<p>7.5cm x 8.5cm 10cm x 12cm 15cm x 20cm 20cm x 50cm</p>	<p>294600 294700 294800 294500</p>

MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

POST-OPERATIVE DRESSINGS

Mepilex® Border Post-Op

Self-adherent post-operative dressing with Flex Technology



SafetaC
TECHNOLOGY

Added benefits:

- Flex Technology provides better conformability
- Increased total contact area in dressing's border

- Low to heavily exuding wounds
- Surgical wounds, cuts and abrasions
- Optimised for post-op use and blood absorption
- Ideal to be used over joints such as hips and knees

6cm x 8cm
9cm x 10cm
10cm x 15cm
10cm x 20cm
10cm x 25cm
10cm x 30cm
10cm x 35cm

496100
496200
496300
496400
496450
496600
496650

Mepore® Film&Pad

Breathable, transparent, absorbent island dressing



- For a wide variety of wounds which have low to moderate exudate levels
- Post-operative wounds

9cm x 20cm
9cm x 25cm
9cm x 35cm

275600
275700
275900

FIBRE / ALGINATE DRESSINGS

Exufiber®

Gelling fibre dressing with Hydrolock® Technology



Exufiber is intended to be used on a wide range of exuding wounds:

- Leg ulcers and diabetic foot ulcers
- Pressure injuries
- Partial thickness burns
- Surgical wounds
- Donor sites
- Malignant wounds

5cm x 5cm
10cm x 10cm
15cm x 15cm
20cm x 30cm
4.5cm x 20cm
1cm x 45cm
(Cavity)
2cm x 45cm
(Cavity)

709900
709901
709903
709904
709906
709908
709909

Melgisorb® Plus

Highly absorbent alginate dressing



- Moderate to heavily exuding wounds
- Pressure injuries
- Venous and arterial leg ulcers
- Diabetic foot ulcers
- Donor sites
- Post-operative wounds
- Dermal lesions and other external wounds inflicted by trauma

5cm x 5cm
10cm x 10cm
10cm x 20cm
3cm x 45cm
(Cavity)

252000
252200
252500
253500

ACTIVE DRESSING

Mesalt®

Sodium chloride impregnated dressing



- For the management of heavily draining and infected wounds in the inflammatory phase
- Deep cavity wounds

5cm x 5cm
7.5cm x 7.5cm
10cm x 10cm
2cm x 1m
(Cavity)

285580
285780
286080
285280

MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

ANTIMICROBIAL DRESSINGS

	<p>Mepilex® Border Ag Antimicrobial absorbent self-adherent foam dressing</p>  	<ul style="list-style-type: none"> Moderate to heavily exuding wounds Infected wounds Leg ulcers Pressure injuries Malignant wounds Partial thickness burns Surgical incisions 	<p>7.5cm x 7.5cm 10cm x 10cm 15cm x 15cm 15cm x 20cm Surgical Sizing 10cm x 20cm 10cm x 25cm 10cm x 30cm</p>	<p>395200 395300 395400 395600 395800 395700 395900</p>
	<p>Mepilex® Border Sacrum Ag Antimicrobial absorbent self-adherent foam dressing for the sacral area</p>  	<ul style="list-style-type: none"> Moderate to heavily exuding wounds for the sacral area where a self-adherent antimicrobial dressing is required 	<p>18cm x 18cm 23cm x 23cm</p>	<p>382000 382400</p>
	<p>Mepilex® Ag Antimicrobial absorbent foam dressing</p>  	<ul style="list-style-type: none"> Low to moderately exuding wounds Infected wounds Leg ulcers and foot ulcers Pressure injuries Partial thickness burns Can be used under compression 	<p>10cm x 10cm 10cm x 20cm 15cm x 15cm 20cm x 20cm 20cm x 50cm</p>	<p>287110 287210 287310 287410 287510</p>
	<p>Mepilex® Heel Ag Antimicrobial absorbent foam dressing for the heel area</p>  	<ul style="list-style-type: none"> Low to moderately exuding wounds Infected wounds Diabetic foot ulcers Pressure injuries Can be used under compression 	<p>13cm x 20cm</p>	<p>388100</p>
	<p>Mepilex® Transfer Ag Antimicrobial compressed foam, exudate transfer dressing</p>  	<ul style="list-style-type: none"> Highly exuding wounds Infected wounds or malignant wounds Skin grafts and donor sites Partial thickness burns Leg ulcers and diabetic foot ulcers Traumatic wounds 	<p>7.5cm x 8.5cm 10cm x 12.5cm 12.5cm x 12.5cm 15cm x 20cm 20cm x 50cm</p>	<p>394000 394100 394700 394800 394500</p>
	<p>Melgisorb® Ag Highly absorbent antimicrobial alginate dressing</p> 	<ul style="list-style-type: none"> Moderate to heavily exuding wounds Post-operative wounds Leg ulcers and diabetic foot ulcers Trauma wounds (dermal lesions, trauma injuries or incisions) Cavity wounds Graft and donor sites Pressure injuries Superficial or partial thickness burns 	<p>5cm x 5cm 10cm x 10cm 15cm x 15cm 20cm x 30cm 3cm x 44cm (Cavity)</p>	<p>256055 256105 256155 256205 256605</p>
	<p>Exufiber® Ag+ Antimicrobial gelling fibre dressing with Hydrolock® technology</p> 	<p>Exufiber Ag+ is intended to be used on a wide range of exuding wounds:</p> <ul style="list-style-type: none"> Leg ulcers and diabetic foot ulcers Pressure injuries Partial thickness burns Surgical wounds Donor sites Malignant wounds 	<p>5cm x 5cm 10cm x 10cm 15cm x 15cm 20cm x 30cm 2cm x 45cm (Cavity)</p>	<p>603421 603425 603423 603424 603420</p>

MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

SUPERABSORBENTS



Mextra® Superabsorbent

Super absorbent dressing with a fluid-repellent backing



- Moderate to very heavily exuding wounds
- Can be used as a primary or secondary dressing where high absorbency is required e.g. together with Mepitel, Mepitel One or Mepilex Transfer
- Minimises the risk of maceration
- Use with Tubifast for light retention

12.5cm x 12.5cm
12.5cm x 17.5cm
12.5cm x 22.5cm
17.5cm x 22.5cm
22.5cm x 27.5cm
22.5cm x 32.5cm
22.5cm x 42.5cm

610000
610100
610200
610300
610400
610500
610600



Mesorb®

Highly absorbent dressing with a fluid-repellent backing



- Moderate to heavily exuding wounds
- Drain sites
- Can be used as a secondary absorbent dressing, where high absorbency is required e.g. together with Mepitel®, Mepitel One or Mepilex® Transfer
- Minimise the risk of maceration
- Use with Tubifast® for light retention

10cm x 13cm
13cm x 15cm
10cm x 23cm
15cm x 23cm
23cm x 25cm
23cm x 30cm

677001
677201
677401
677601
677701
677801

FILM DRESSINGS



Mepitel® Film

Breathable, transparent film dressing with Safetac® silicone adhesive.



- Superficial skin injuries and burns
- Protects fragile and sensitive skin
- Protective cover for open surgical wounds such as abdomen
- To fixate primary dressings
- Can be used in combination with gels and ointments

6cm x 7cm
10cm x 12cm
10cm x 25cm
15cm x 20cm

296170
296270
296470
296670



Mepore® Film

Breathable, transparent film dressing

- Designed as a secondary fixation dressing and as a primary dressing for a wide range of clean wounds

6cm x 7cm
10cm x 12cm
10cm x 25cm
15cm x 20cm

270600
271500
272500
273000

ISLAND DRESSINGS



Mepore®

Absorbent island dressing



- For a wide variety of superficial wounds which have low to moderate levels of exudate
- Post-operative wounds
- Cuts and grazes

6cm x 7cm
9cm x 10cm
9cm x 15cm
9cm x 20cm
9cm x 25cm
9cm x 30cm
9cm x 35cm

670800
670900
671000
671100
671200
671300
671400



Mepore® Pro

Shower-proof, absorbent island dressing



- For a wide variety of superficial wounds which have low to moderate levels of exudate where a shower-proof dressing is required
- Post-operative wounds
- Cuts and grazes

6cm x 7cm
9cm x 10cm
9cm x 15cm
9cm x 20cm
9cm x 25cm
9cm x 30cm

670820
670920
671020
671120
671220
671320



Mepore® Film&Pad

Transparent, absorbent island dressing



- For a wide variety of wounds which have low to moderate exudate levels
- Cuts and abrasions
- Post-operative wounds
- Minor burns

4cm x 5cm
5cm x 7cm
(Oval)
5cm x 7cm
(Square)
9cm x 10cm
9cm x 15cm
9cm x 20cm
9cm x 25cm
9cm x 35cm

275100
275200
275300
275400
275500
275600
275700
275900

MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

BANDAGES & DERMATOLOGY



Tubifast® 2-WAY WITH STRETCH TECHNOLOGY

Lightweight tubular dressing retention bandage

- Lightweight tubular dressing retention bandage that is ideal for any part of the body
- Can be used for patch wrapping and as an undercast stockinette

10m roll
RED
 3.5cm x 10m
GREEN
 5cm x 10m
BLUE
 7.5cm x 10m
YELLOW
 10.75cm x 10m
PURPLE
 20cm x 10m

2434
 2436
 2438
 2440
 2444



Tubifast® Garments

Ready-to-wear garments

- Ready-made garments, for dressing retention and wet or dry wrapping, are available in full-sleeved vests, tights, leggings, socks and gloves
- For children of six months to 14 years of age
- With no tapes, pins or ties, Tubifast Garments are quick and easy to apply, and their two-way stretch construction ensures complete freedom of movement, so aiding compliance. Comfortable to wear under nightwear and ordinary clothes

Vests
 6-24 months 992007
 2-5 yr 992008
 5-8 yr 992009
 8-11 yr 992010
 11-14 yr 992011
Tights
 6-24 months 992012
Leggings
 2-5 yr 992013
 5-8 yr 992014
 8-11 yr 992015
 11-14 yr 992016
Socks
 1 size: 2-14yr 992017
Gloves
 Small child 5923



Tubigrip®

Multi-purpose, elasticated tubular support bandage for treatment of strains and sprains, soft tissue injuries, general oedema, post-burn scarring and rib cage injuries. Can also be used for pressure dressings and arm fixation

Assortment (All 10m)

- A - infant legs and arms
- B - small wrists
- C - medium wrists, small elbows
- D - large wrists, medium elbows, small/medium ankles, medium knees
- E - large elbows, medium/large ankles, medium knees
- F - large knees, small thighs
- G - x large knees, medium thighs
- J - large thighs, small trunks
- K - small/medium trunks
- L - large trunks
- M - X large trunks

Natural	Beige
1435	-
1436	1449
1443	1450
1437	1451
1434	1448
1438	1452
1439	1453
1440	-
1441	-
1442	-
1444	-



Tubinette®

Lightweight tubular dressing retention bandage

- Tubinette holds dressings securely and comfortably in place, without constriction or compression. The tubular construction makes application quick and easy

All 20m
 01 Fingers&toes 2416
 12 Bulky fingers 2417

MÖLNLYCKE - WOUND CARE PRODUCTS

PRODUCT DESCRIPTION

AREAS OF USE

SIZE

CODE

BANDAGES & DERMATOLOGY



Setopress®
High compression bandage

- Used in the treatment of venous leg ulcers and associated conditions requiring graduated compression of the leg.

10cm x 3.5m

3505



Elset® S / Elset®
Light compression bandages

- Used to provide compression for the leg following vascular surgery and for general surgical and orthopaedic use
- Specifically recommended for above and below knee stump bandaging

Elset 'S'

10cm x 12m

15cm x 12m

Elset

10cm x 6m

10cm x 8m

993486

993487

3482

3489

FIXATION



Mepitac®
Polyurethane tape with **Safetac®** technology

Safetac
TECHNOLOGY

- Fixation of medical devices such as drains, tubes and probes
- Skin protection under devices

2cm x 3m

4cm x 1.5m

298300

298400



Mefix®
Self-adherent fabric fixation tape

- For fixing a wide variety of dressings, swabs, catheters and cannulae

10m roll

x 2.5cm

x 5cm

x 10cm

x 15cm

x 20cm

x 30cm

310299

310599

311099

311599

312000

313000

SCAR THERAPY



Mepiform®
Scar dressing with **Safetac** silicone adhesive

Safetac
TECHNOLOGY

- Keloid scars
- Hypertrophic scars

4cm x 30cm

5cm x 7.5cm

10cm x 18cm

293100

293200

293400

SKIN THERAPY



Epaderm® Cream
An emollient for dry skin conditions

- A 2-in-1 emollient - it can be used as a skin emollient as well as a skin cleanser
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages, including babies

25g

50g

500g

99400822

99400851

99400850



Epaderm® Ointment
An emollient for dry skin conditions

- Epaderm Ointment is a 3-in-1 emollient
- Used for eczema, psoriasis and all other dry skin conditions
- Suitable for all ages, including babies
- Can be used on the skin, as a bath additive or a skin cleanser

25g

125g

500g

1kg

400861

99400846

99400847

99400818

MÖLNLYCKE - WOUND CARE PRODUCTS




PRODUCT DESCRIPTION

AREAS OF USE


SIZE

CODE

OTHER DRESSINGS

	<p>Mepore® IV Transparent, breathable IV film dressing</p>	<ul style="list-style-type: none"> Mepore IV is intended for fixation of IV catheters Mepore IV is <u>not</u> intended for fixation of life supporting devices 	<p>5.5cm x 5cm 8cm x 9cm 10cm x 11cm</p>	<p>274000 274200 274400</p>
	<p>Lyof foam® Max T Polyurethane foam dressing for use around tracheostomy tubes</p> 	<ul style="list-style-type: none"> A soft tracheostomy dressing that absorbs exudate / secretion and provides protection of the skin around the tracheostomy tube 	<p>9cm x 9cm</p>	<p>603207</p>

TURNING & POSITIONING SYSTEMS

	<p>Mölnlycke® Z-Flex® Heel Boot To maintain offloading of heels overtime. The combination of positive air displacement and fluidised positioning helps offload the heel and spread the pressure evenly over the Achilles tendon.</p>	<ul style="list-style-type: none"> Fully offloads the heel Positions feet in neutral alignment Accommodates compression devices Quick access for heel inspection Thermal regulating material One size fits most 	<p>One size fits most 2 boots per box 8 boots per box</p>	<p>1400122 1400123</p>
	<p>Mölnlycke® Z-Flo® Fluidized Positioner To assist with patient positioning and offloading of at risk anatomical sites.</p>	<ul style="list-style-type: none"> Redistributes pressure or offloads Keeps patients in position Mouldable to protect vulnerable areas 	<p>With Tab 41cm x 76cm Without Tab 64cm x 91cm 29cm x 50cm (Occipital)</p>	<p>1401003 1401005 1401007</p>
	<p>Mölnlycke® Tortoise® Turning & Positioning System To assist with turning and repositioning of patients.</p>	<ul style="list-style-type: none"> Provides continuous pressure redistribution Supports a wide degree of turns to increase efficiency Can be used to create any degree of turn required by the patient's condition Reduces the risk of caregiver injury Reduces the risk of the patient falling or sliding down in bed Approved for lateral transfer and MRI to minimise lifting 	<p>Standard Bariatric</p>	<p>1400800 1400801</p>

NB: Please refer to Instructions of Use prior to usage.

Find out more at www.molnlycke.co.nz

The Mölnlycke, Mepilex, Safetac, Mepitel, Melgisorb, Mepitac, Mefix, Mepiform, Mesalt, Exufiber, Mepore, Mesorb, Mextra, Tubifast, Tubinette, Epaderm, Hydrolock, Lyof foam, Tubigrip, TSSB, Elset, Setopress, 2-Way Stretch, Z-Flo, Tortoise names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Z-Flex is a trademark used by Mölnlycke Health Care Group of Companies. ©2020 Mölnlycke Health Care. All rights reserved. www.molnlycke.com.au. Mölnlycke Health Care, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231. www.molnlycke.co.nz. AUWC_11122017_V3.ANZ